

PRAXIFLOW

FABRYCZNIE NAPEŁNIONE STRZYKAWKI DO PŁUKANIA I ZABEZPIECZENIA URZĄDZEŃ
DOSTĘPU NACZYNIOWEGO

Każdy zestaw PraxiFlow™ składa się z jednej strzykawki Praxiject™ 0.9% NaCl do płukania i jednej strzykawki CitraFlow™ cytrynianu sodu do zabezpieczenia urządzeń dostępu naczyniowego. Wszystkie udowodnione korzyści z produktów Praxiject™ i CitraFlow™ w fabrycznie napełnionych strzykawkach, w jednym opakowaniu wygodnym w stosowaniu.

- o Opakowanie łatwe do otwierania. Przeznaczone do stosowania w polu jałowym.
- o Łatwe do stosowania fabrycznie napełnione strzykawki. Oszczędność czasu i mniej manipulacji.
- o Każde opakowanie zawiera jedną strzykawkę VAD do płukania i jedną strzykawkę VAD do zabezpieczenia
- o Mniejszy wpływ na środowisko. Mniejsze zużycie opakowań niż w przypadku pojedynczo sprzedawanych produktów.
- o Produktu dostępny w strzykawkach o pojemności 10 ml do płukania oraz różnych stężeniach cytrynianu sodu do zabezpieczenia.

Catalog #	Description	Quantity/case
3215	1 Praxiject™ 0.9% NaCl syringe. (10ml in 10cc syringe) 1 CitraFlow™ 46.7% Sodium Citrate syringe. (3ml in 5cc syringe)	75 packages / cs
3225	1 Praxiject™ 0.9% NaCl syringe. (10ml in 10cc syringe) 1 CitraFlow™ 30% Sodium Citrate syringe. (3ml in 5cc syringe)	75 packages / cs
3255	1 Praxiject™ 0.9% NaCl syringe. (10ml in a 10cc syringe) 1 CitraFlow™ 4% Sodium Citrate syringe. (3ml in 5cc syringe)	75 packages / cs